

ABRF 2026 ANNUAL MEETING



INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION


PITTSBURGH, PA • MARCH 28-31, 2026



UNIVERSITY of MARYLAND
MARLENE AND STEWART GREENEBAUM
COMPREHENSIVE CANCER CENTER



ISAC RECOGNIZED
Shared Resource Laboratory



No Cells Left Behind: Wrangling the Full Biological Discovery Potential out of Spectral Flow Cytometry (SFC) datasets

by David Rach

UMGCC Flow Cytometry Shared Resource

March 30, 2026




INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026



Spectral Flow Cytometry (SFC) is awesome!



- Capacity to resolve spectrally similar fluorophores: greater number markers per panel than Conventional Flow Cytometry (CFC)
- Faster acquisition speed: typically far more cells per .fcs file than Mass Cytometry (MC)
- Analytical breadth and depth permits characterization rare cell populations (and potentially millions of other non-target cells)
- Invaluable for human clinical samples with limited biospecimen




INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026




Complex Datasets, Combinatorially Exponential Clusters

- If you are interested in a single population of cells, analysis process remains essentially unchanged.
- If your goal is to characterize all cell populations that may be present, the breadth and depth complicate the analysis.



INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026




Complex Datasets, Combinatorially Exponential Clusters

- If you are interested in a single population of cells, analysis process remains essentially unchanged.
- If your goal is to characterize all cell populations that may be present, the breadth and depth complicate the analysis.


If we step back, and consider for a given marker:

- We split cells positive and negative populations
- Repeat the splitting process for every marker
- Group cells that have shared expression of every single marker
- Theoretical upper number of clusters would be 2^n



INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026




Complex Datasets, Combinatorially Exponential Clusters

- If you are interested in a single population of cells, analysis process remains essentially unchanged.
- If your goal is to characterize all cell populations that may be present, the breadth and depth complicate the analysis.

If we step back, and consider for a given marker:

- We split cells positive and negative populations
- Repeat the splitting process for every marker
- Group cells that have shared expression of every single marker
- Theoretical upper number of clusters would be 2^n



- **9-color** CFC panel: **512** clusters at most
- **40-color** SFC panel: **1.09e+12** clusters at most




INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION


PITTSBURGH, PA • MARCH 28-31, 2026




A crossroads deal with the Spectral Signature Devil



Unmixing Controls
(Single Color and
Unstained)





Single cell from
full-stained sample



What fluorophores
are present, and
roughly how much


- Consequently, there is a need for semi-supervised and unsupervised algorithmic analytical approaches
- Immense potential for biological discovery, but at the cost of a major trade-off
- Quality of our final data depends on our ability to infer the fluorophores present on individual fully-stained cells
- How closely the fluorescent signatures of our unmixing controls (single-color and unstained) match the spectral signatures as present on the full-stained cells is critical to this final data quality.

“We discovered this new dendritic cell cluster....”




INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026



More unmixing errors, more downstream batch effects




- Despite the critical importance unmixing controls play, few tools permit us to determine quality of the individual control before unmixing occurs
- Existing quality control algorithms (PeacoQC, FlowAI, etc.) primarily focus on detecting laser and fluidic issues in already acquired .fcs files
- Unmixing errors can be introduced by tandem degradation, additional autofluorescences, undiagnosed instrumental issues, etc.
- Unmixing errors contribute to batch effects, and require complicated troubleshooting to be able to proceed with the downstream analysis



INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026





Open-source tool sets specifically for Spectral Flow Cytometry's quirks



Spectral Signature
Characterization and Isolation



Automated Instrument
Quality-Control Dashboards




Semi-supervised analysis via
incorporation of manually-
annotated gating information



INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026






Spectral Signature Characterization and Clustering



- For individual cells, it characterizes the cell's normalized fluorescent signature
- Cells that share normalized spectral signature peaks at individual detectors are then clustered
- **Fluorophore vials in good condition**: we retrieve few normalized signatures from the raw .fcs file for each single-color unmixing control
- When the **tandem fluorophore decouples or degrades**, noticeable increase in variant signatures, flag the issue
- Setting up automation to retrieve and process reference controls for individual users experiments during off-hours



INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026











Automated isolation of autofluorescence signatures



- Adapt this approach to isolate out
- normalized autofluorescence signatures
- from unstained controls


- Create .fcs files with purified autofluorescence signatures, to use for unmixing with commercial software.

- Isolated cells share same normalized signature, only differing in brightness, unlike the current Autofluorescence Explorer or gating a dimensionality-visualized island style approaches to autofluorescence isolation






INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026





Residual Autofluorescence peaks in Single-Color Signatures




For cell controls with a **moderate density antigen** on **bright fluorophores**, we see that any residual autofluorescence contributes minimally to the overall isolated single-color signature when gated according to best practices






INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026





Residual Autofluorescence peaks in Single-Color Signatures




For cell controls with a **moderate density antigen** on **dim fluorophores**, we see substantial residual autofluorescence contributions to the isolated single-color normalized signatures even when gated according to best practices



INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026







Residual Autofluorescence peaks in Single-Color Signatures



**PE signature with
residual AF**



**PE signature without
residual AF**




This residual autofluorescence contribution results in unmixing errors, and the panel unmixing significantly improved when switched out for bead controls




INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026






Instrument Quality-Control Dashboards




- Prior work assumes that the spectral flow cytometers are stably maintaining their quality control (QC) targets
- At our core: Daily acquire 5000 QC beads as .fcs files both before and after morning QC
- Automatically processed for MFI, Gain, and %rCV and displayed as Levey-Jennings plots
- Data for all instruments is deployed as interactive plots on a Quarto website hosted on GitHub





INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026





Early detection of Instrumental Issues




- Permitted identification of laser issues a month before standard QC would have failed, allowing earlier scheduling of preventive maintenance visits
- By spotting the spike in %rCV values for last detector of each laser, we were able to link unmixing errors we were seeing to the Lot 2006 QC beads being bleach-sensitive



INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026






When automated QC goes awry




- Bleach stripped APC-Fire 810 like fluorophore from the QC beads
- MFI for the detectors where fluorophore signal are below the target value, gain is increased by automated QC
- Major unmixing issues result, because the instrument settings didn't need the gain increase
- Retrospective analysis suggest issue has been present for other bead lots, but was never flagged.



INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026



Coereba




Semi-supervised analysis by annotating in manual gates




- Comprehensive automated manual-gating information is appended to individual cells in an .fcs file
- Can be passed to unsupervised algorithms for quality control, normalization, clustering, dimensionality visualization
- Can retrieve underlying supervised gating information afterwards

Marker	MAITs	NKTs	Vδ2 CD4	Vδ2 CD8	Vδ2 DN	Vδ2 T17	Vδ2 PD1
CD7	+++	+++	+++	+++	+++	+++	+++
CD3	+++	+++	+++	+++	+++	+++	+++
Vδ2	-	-	+++	+++	+++	+++	+++
Vα7.2	+++	-	-	-	-	-	-
Vα24α18	-	+++	-	-	-	-	-
CD8	++	-	-	+++	-	-	++
CD4	-	+++	+++	-	-	-	-
CD26	+++	++	++	-	-/+	+++	-
CCR6	+++	-	-	-	-	+++	-
CD161	+++	++	++	++	+	+++	-
NKG2D	+++	-	-	+++	+++	+++	-
CD56	-	-	-	-/+	-	-	-/+
CD16	-	-	-	+	-	-	-/+
CD62L	-	-	-	-/+	-	-	-
CD27	+++	+++	+++	+++	+++	+++	+++
CD45RA	+++	-	-	+++	++	+++	+++
CD127	+++	+++	+++	+	+++	+++	-
CD25	+	+++	+++	-	-	+++	-
CCR4	-	++	++	-	-	+	-
PD1	+	++	++	++	++	+	+++
CD38	-	+	-/+	+	+	-	-




INNOVATING AT THE INTERSECTION OF
 SCIENCE, TECHNOLOGY, AND COLLABORATION
 PITTSBURGH, PA • MARCH 28-31, 2026



Current Status

- All packages are available via our GitHub
- Fixing bugs and updating documentation
- Submitting Bioconductor this year




INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026



Current Status

- All packages are available via our GitHub
- Fixing Bugs and updating Documentation
- Submitting Bioconductor this year



```
— R CMD build —
✓
checking for file '/tmp/Rtmpo0G4tw/remotes13712b1288d8/DillonHamill-CytoExploreR-71a8da6/D
DESCRIPTION'
...
- preparing 'CytoExploreR':
✓ checking DESCRIPTION meta-information ...
- checking for LF line-endings in source and make files and shell scripts
- checking for empty or unneeded directories
  Omitted 'LazyData' from DESCRIPTION
- building 'CytoExploreR_1.1.0.tar.gz' (475ms)

Installing package into '/home/david/R/x86_64-pc-linux-gnu-library/4.5'
(as 'lib' is unspecified)
ERROR: dependency 'EmbedSOM' is not available for package 'CytoExploreR'
Perhaps try a variation of:
install.packages('EmbedSOM')
* removing '/home/david/R/x86_64-pc-linux-gnu-library/4.5/CytoExploreR'
Warning message:
In i.p(...) :
  installation of package '/tmp/Rtmpo0G4tw/file137111e5ea27/CytoExploreR_1.1.0.tar.gz' had
non-zero exit status
>
```




- What good is a tool that the community can't use?



INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026






UNIVERSITY of MARYLAND
MARLENE AND STEWART GREENEBAUM
COMPREHENSIVE CANCER CENTER

Flow Cytometry
Shared Resource


ISAC RECOGNIZED
Shared Resource Laboratory

CYTOMETRY IN R

A COURSE FOR BEGINNERS




- Free weekly Cytometry in R course, offered both in-person and online via YouTube livestreams.
- Aimed at those with prior cytometry experience, but no-to-little coding experience
- Started in February, will be exploring from beginner to advanced topics for at least 30 weeks
- Around 2000 interest forms completed, between 300-500 weekly participants



INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026



Acknowledgements



ISAC RECOGNIZED
Shared Resource Laboratory



UMGCCC FCSR

Dr. Xiaoxuan Fan
Dr. Natarajan Ayithan
Mikayla Trainor

- University of Maryland School of Medicine

- Dr. Cristiana Cairo
- Dr. Kirsten E. Lyke

PIN:


Daniela
Franco

U01 HD092308


This publication was supported by funds through Maryland Department of Health's Cigarette Restitution Fund Program – CH- 649-CRF and the National Cancer Institute – Cancer Center Support Grant (CCSG) – P30CA134274.d Service

- Molecular Microbiology Immunology Graduate Program

- Dr. Eileen Barry
- Bess Tracey



Cytometry
Discord



r/flowcytometry

Bioconductor Zulip



<https://github.com/DavidRach>




@davidrach.bsky.social




INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026




Questions?



UMGCCC FCSR
InstrumentQC dashboard



Luciernaga




Website - CytometryInR



InstrumentQC How-To



Coereba



YouTube @CytometryInR



INNOVATING AT THE INTERSECTION OF
SCIENCE, TECHNOLOGY, AND COLLABORATION

PITTSBURGH, PA • MARCH 28-31, 2026

